

UC SANTA BARBARA

# *THE* **Current**

March 7, 2011

Gail Gallessich

## **New Instrument Keeps An Electronic 'Eye' on Nanoparticles**

Precision measurement in the world of nanoparticles has now become a possibility, thanks to scientists at UC Santa Barbara.

The UCSB research team has developed a new instrument capable of detecting individual nanoparticles with diameters as small as a few tens of nanometers. The study will be published on line this week by Nature Nanotechnology, and appear in the April print issue of the journal.

"This device opens up a wide range of potential applications in nanoparticle analysis," said Jean-Luc Fraikin, the lead author on the study. "Applications in water analysis, pharmaceutical development, and other biomedical areas are likely to be developed using this new technology." The instrument was developed in the lab of Andrew Cleland, professor of physics at UCSB, in collaboration with the group of Erkki Ruoslahti, Distinguished Professor, Sanford-Burnham Medical Research Institute at UCSB.

Fraikin is presently a postdoctoral associate in the Marth Lab at the Sanford-Burnham Medical Research Institute's Center for Nanomedicine, and in the Soh Lab in the Department of Mechanical Engineering at UC Santa Barbara.

The device detects the tiny particles, suspended in fluid, as they flow one by one through the instrument at rates estimated to be as high as half a million particles

per second. Fraikin compares the device to a nanoscale turnstile, which can count -- and measure -- particles as they pass individually through the electronic "eye" of the instrument.

The instrument measures the volume of each nanoparticle, allowing for very rapid and precise size analysis of complex mixtures. Additionally, the researchers showed that the instrument could detect bacterial virus particles, both in saline solution as well as in mouse blood plasma.

In this study, the researchers further discovered a surprisingly high concentration of nanoparticles present in the native blood plasma. These particles exhibited an intriguing size distribution, with particle concentration increasing as the diameter fell to an order of 30 to 40 nanometers, an as-yet unexplained result.

---

[\[RETURN TO TOP\]](#)

† Top photo: Andrew Cleland (left) and Jean-Luc Fraikin

Credit: George Foulsham, Office of Public Affairs, UCSB

†† Middle photo: Optical microscope image of the microfluidic channel (light pattern) and sensing electrode (gold) of the analyzer.

Nanoparticles suspended in a fluid flow through the channel and are detected individually as they pass through the sensing volume.

Credit: J.L. Fraikin and A.N. Cleland, UCSB

††† Bottom photo: Scanning electron microscope image of polystyrene nanoparticles. These particles are roughly 100 nanometers in diameter, and are readily detected with the analyzer. The uniform size distribution is used to calibrate the instrument.

Credit:

J.L. Fraikin and A.N. Cleland, UCSB

Related Links

## **About UC Santa Barbara**

The University of California, Santa Barbara is a leading research institution that also provides a comprehensive liberal arts learning experience. Our academic community of faculty, students, and staff is characterized by a culture of interdisciplinary collaboration that is responsive to the needs of our multicultural and global society. All of this takes place within a living and learning environment like no other, as we draw inspiration from the beauty and resources of our extraordinary location at the edge of the Pacific Ocean.